



EXCLUSIVE
SensorDisc
TECHNOLOGY

■ **LONG LIFE
BATTERY**

RoHS

NEWLINE[®]

Electronic Talking Bathroom Scale

INSTRUCTION

For the purpose of proper using & maintaining the scale, please read this instruction completely before operating.

NEWLINE®

Electronic Talking Bathroom Scale

Your new electronic personal scale has been designed to accurately indicate your weight gain or loss over a period of time and should give many years of service under normal use .

SPECIFICATIONS

- Equipped with a high precision strain gauge sensor system
- Capacity: 150kgs 180 kgs
- Division: 0.1kg
- Voice indicate the weight data
- Foot-tap switch
- Automatic zero resetting / switch off
- Unit switch: kg / lb / st
- Low battery / over-load indicator

POWER SUPPLY

- Battery 2xCR2032 1x9V 2xAAA
- Insert the AAA battery before the first use.
- When the LCD indicate "Lo", open the battery door and replace the battery with a new one.



OPERATION

STEP1: Put the scale on a hard & flat surface (avoid carpet or soft surface)

STEP2: Switch on the scale by slightly tapping the middle front platform, then wait a few seconds, the following will appear :





STEP3: Gently step on the scale. Put your feet evenly on the scale without moving and wait while the scale computes your weight. Your actual weight will be displayed on the screen.



VOICE FUNCTION

1. Stand on the center of scale stillly, then the LCD will show the data and at the same time, the scale will voice the data.
2. You can adjust the switch to control the voice volume.



AUTOMATIC SWITCH OFF

When you step off the scale, it will automatically switch off, (if no weighing is carried out while the zero is being displayed, the scale switch itself off automatically after approximate 6 seconds.)

OVER-LOAD INDICATOR

When the scale is over loaded (over 150kg 180kg) "Err" or "OL" will appear on the screen:



ADVICE FOR USE AND CARE

- Placing your scale on a hard & flat surface will ensure the greatest accuracy and repeatability.
- Always weigh yourself on the same scale placed on the same floor surface. Uneven floors may affect readings.
- Weigh yourself without clothing and footwear; before meals and always at the same time of day. The best time to weigh yourself is in the morning.
- Your scale rounds up or down to the nearest 100g, if you weigh yourself twice and get two different readings; your weight lies between the two.
- Clean your scale with a damp cloth but do not allow water to enter the scale. Do not use chemical cleaning agents.
- Do not allow your scale to become saturated with water, this can damage the electronics.
- Do not store the scale in an upright position when the scale is not in use.
- Remove the batteries if the scale is not being used for a long time . Treat your scale with care. It is a precision instrument. Do not drop it or jump on it.
- Do not use it for legal trade.